


Product Datasheet

HLA-DRB5 antibody (orb539833)

Description	Rabbit polyclonal antibody to HLA-DRB5
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human HLA-DRB5(NP_002116.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	37kDa
Uniprot ID	Q30154
Dilution Range	WB: 1:1000-1:3000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...