



Product Datasheet

HMGA1 antibody (orb539838)


Description	Rabbit polyclonal antibody to HMGA1
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of human HMGA1(NP_665908.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	16kDa
Uniprot ID	P17096
Dilution Range	WB: 1:1000-1:2000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse kidney tissue.



Immunohistochemical staining of mouse brain tissue.



Western blot analysis of extracts of rat kidney.