


Product Datasheet

HYAL3 antibody (orb539900)

Description	Rabbit polyclonal antibody to HYAL3
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant fusion protein of human HYAL3 (NP_003540.2).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	47kDa
Uniprot ID	043820
Entrez	8372
Dilution Range	WB: 1:500-1:2000, IF: 1:50-1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...