



Product Datasheet

IDH1 antibody (orb539911)


Description	Rabbit polyclonal antibody to IDH1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IP, WB
Immunogen	Recombinant fusion protein of human IDH1(NP_005887.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	46KDa
Uniprot ID	O75874
Dilution Range	WB: 1:500-1:2000, IP: 1:50-1:100
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of L929 cell...



Immunofluorescence analysis of C6 cells



Western blot analysis of extracts of var...