


## Product Datasheet

### IDI1 antibody (orb539915)

|                            |                                                        |
|----------------------------|--------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to IDI1                     |
| <b>Reactivity</b>          | Human, Mouse                                           |
| <b>Conjugation</b>         | Unconjugated                                           |
| <b>Tested Applications</b> | WB                                                     |
| <b>Immunogen</b>           | Recombinant fusion protein of human IDI1(NP_004499.2). |
| <b>Preservatives</b>       | PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.       |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.            |
| <b>Note</b>                | For research use only                                  |
| <b>Isotype</b>             | IgG                                                    |
| <b>Clonality</b>           | Polyclonal                                             |
| <b>MW</b>                  | 26kDa                                                  |
| <b>Uniprot ID</b>          | <b>Q13907</b>                                          |
| <b>Dilution Range</b>      | WB: 1:500-1:2000                                       |
| <b>Expiration Date</b>     | 12 months from date of receipt.                        |



Western blot analysis of extracts of var...