



Product Datasheet

IGF2BP2 antibody (orb539940)


Description	Rabbit polyclonal antibody to IGF2BP2
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human IGF2BP2/IMP2(NP_006539.3).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	70kDa
Uniprot ID	Q9Y6M1
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of rat testis...



Immunofluorescence analysis of U-2 OS cells...



Western blot analysis of extracts of various cell lines...