



Product Datasheet

IGF2BP3 antibody (orb539942)

Description	Rabbit polyclonal antibody to IGF2BP3
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, IP, WB
Immunogen	Recombinant fusion protein of human IGF2BP3(NP_006538.2).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	64KDa
Uniprot ID	000425
Dilution Range	WB: 1:500-1:1000, IF: 1:50-1:100, IP: 1:20-1:50
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of U2OS cell...



Western blot analysis of extracts of var...