

Product Datasheet

DSG3 antibody (orb53995)

Description	DSG3 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Desmoglein-3 protein (50-615AA)
Target	DSG3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P32926
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse brain tis...