



## Product Datasheet

### IL18 antibody (orb539965)

|                            |                                                 |
|----------------------------|-------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to IL18              |
| <b>Reactivity</b>          | Human, Mouse, Rat                               |
| <b>Conjugation</b>         | Unconjugated                                    |
| <b>Tested Applications</b> | ICC, IF, IP, WB                                 |
| <b>Immunogen</b>           | A synthetic peptide of human IL18(NP_001553.1). |
| <b>Preservatives</b>       | PBS with 0.05% proclin300,50% glycerol,pH7.3.   |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.     |
| <b>Note</b>                | For research use only                           |
| <b>Isotype</b>             | IgG                                             |
| <b>Clonality</b>           | Polyclonal                                      |
| <b>MW</b>                  | 18KDa/22KDa                                     |
| <b>Uniprot ID</b>          | <b>Q14116</b>                                   |
| <b>Dilution Range</b>      | WB: 1:500-1:2000, IF: 1:50-1:200                |
| <b>Expiration Date</b>     | 12 months from date of receipt.                 |



Immunofluorescence analysis of U2OS cell...



Western blot analysis of extracts of var...