



Product Datasheet

IL6 antibody (orb539985)


Description	Rabbit polyclonal antibody to IL6
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human IL6 (NP_000591.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	23kDa
Uniprot ID	P05231
Entrez	3569
Dilution Range	WB: 1:200-1:500, IHC-P: 1:50-1:200, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human pl...



Immunohistochemical staining of human to...



Western blot analysis of extracts of var...