


Product Datasheet


IQGAP1 antibody (orb540010)

Description	Rabbit polyclonal antibody to IQGAP1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human IQGAP1(NP_003861.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	189kDa
Uniprot ID	P46940
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of C6 cells

...



Western blot analysis of extracts of var...