



## Product Datasheet

### CD11c antibody (orb540031)


|                            |                                                         |
|----------------------------|---------------------------------------------------------|
| <b>Description</b>         | Rabbit Polyclonal antibody to CD11c                     |
| <b>Reactivity</b>          | Mouse, Rat                                              |
| <b>Conjugation</b>         | Unconjugated                                            |
| <b>Tested Applications</b> | ICC, IF                                                 |
| <b>Immunogen</b>           | Recombinant fusion protein of human ITGAX(NP_000878.2). |
| <b>Preservatives</b>       | PBS with 0.02% sodium azide,50% glycerol,pH7.3.         |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.             |
| <b>Note</b>                | For research use only                                   |
| <b>Isotype</b>             | IgG                                                     |
| <b>Clonality</b>           | Polyclonal                                              |
| <b>MW</b>                  | 127kDa                                                  |
| <b>Uniprot ID</b>          | <b>P20702</b>                                           |
| <b>Dilution Range</b>      | WB: 1:500-1:2000, IHC-P: 1:50-1:200                     |
| <b>Expiration Date</b>     | 12 months from date of receipt.                         |



Immunohistochemical staining of human lu...



Immunohistochemical staining of human to...



Immunohistochemical staining of rat lung...