



Product Datasheet

ITPR1 antibody (orb540044)

Description	Rabbit polyclonal antibody to ITPR1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, IP, WB
Immunogen	A synthetic peptide of human IP3R(XP_011531985.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	320kDa
Uniprot ID	Q14643
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200, IP: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of U-2 OS cells...



Western blot analysis of extracts of various tissues...