



Product Datasheet

JAG2 antibody (orb540052)

Description	Rabbit polyclonal antibody to JAG2
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human JAG2(NP_002217.3).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	150kDa
Uniprot ID	Q9Y219
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.




Immunofluorescence analysis of L929 cell...



Immunofluorescence analysis of C6 cells

...



Western blot analysis of extracts of var...