

Product Datasheet

NQO1 antibody (Biotin) (orb54006)

Description	NQO1 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human NAD(P)H dehydrogenase [quinone] 1 protein (55-155AA)
Target	NQO1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P15559
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.