



Product Datasheet

MYST1 antibody (orb540069)

Description	Rabbit polyclonal antibody to MYST1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human KAT8/MYST1/MOF(NP_892003.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	53KDa
Uniprot ID	Q9H7Z6
Dilution Range	WB: 1:200-1:500, IF: 1:100-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of U-2 OS cells



Western blot analysis of extracts of various cell lines