



Product Datasheet

KDM1A antibody (orb540100)


Description	Rabbit polyclonal antibody to KDM1A
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Recombinant fusion protein of human KDM1 (NP_055828.2).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	110kDa
Uniprot ID	O60341
Entrez	23028
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of L929 cell...



Immunofluorescence analysis of HeLa cell...



Western blot analysis of extracts of var...