



Product Datasheet

KIAA1429 antibody (orb540114)


Description	Rabbit polyclonal antibody to KIAA1429
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human KIAA1429(NP_056311.2).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	210KDa
Uniprot ID	Q69YN4
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse brain tissue.



Immunohistochemical staining of human placental tissue.



Immunohistochemical staining of human testis tissue.