



Product Datasheet

KIF1C antibody (orb540117)


Description	Rabbit polyclonal antibody to KIF1C
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human KIF1C(NP_006603.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	165kDa
Uniprot ID	O43896
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse heart tissue.



Immunohistochemical staining of rat heart tissue.



Western blot analysis of extracts of various cell lines.