



Product Datasheet

KIF5B antibody (orb540124)


Description	Rabbit polyclonal antibody to KIF5B
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human KIF5B(NP_004512.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	110KDa
Uniprot ID	P33176
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse pa...



Immunohistochemical staining of mouse sp...



Immunohistochemical staining of mouse te...