



Product Datasheet

KLF5 antibody (orb540141)

Description	Rabbit polyclonal antibody to KLF5
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human KLF5(NP_001721.2).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	55kDa
Uniprot ID	Q13887
Dilution Range	WB: 1:200-1:1000, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of A375 cell...



Western blot analysis of extracts of var...