



Product Datasheet

KLRC4 antibody (orb540158)


Description	Rabbit polyclonal antibody to KLRC4
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human KLRC4(NP_038459.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	14kDa
Uniprot ID	O43908
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of mouse spleen tissue



Immunofluorescence analysis of Jurkat cells



Western blot analysis of extracts of various tissues