



Product Datasheet

KRT17 antibody (orb540175)


Description	Rabbit polyclonal antibody to KRT17
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	A synthetic peptide of human Cytokeratin 17 (KRT17)(NP_000413.1).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	48KDa
Uniprot ID	Q04695
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HeLa cell...



Immunohistochemical staining of rat kidn...



Western blot analysis of extracts of var...