


Product Datasheet

LEFTY1 antibody (orb540207)

Description	Rabbit polyclonal antibody to LEFTY1
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human LEFTY1(NP_066277.1).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	39kDa
Uniprot ID	075610
Dilution Range	WB: 1:200-1:2000
Expiration Date	12 months from date of receipt.



Western blot
analysis of
extracts of
var...