

Product Datasheet

mpt64 antibody (FITC) (orb54021)

Description	mpt64 antibody (FITC)
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	FITC
Immunogen	Recombinant Mycobacterium tuberculosis Immunogenic protein MPT64 protein (24- 228AA)
Target	mpt64
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P9WIN9
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.