



Product Datasheet

Galectin 3 antibody (orb540217)


Description	Rabbit Polyclonal antibody to Galectin 3
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human Galectin 3/Galectin 3/LGALS3(NP_002297.2).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	28KDa
Uniprot ID	P17931
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of U-2 OS cells



Immunofluorescence analysis of L-929 cells



Immunofluorescence analysis of C6 cells