


Product Datasheet

MAP3K13 antibody (orb540316)

Description	Rabbit polyclonal antibody to MAP3K13
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human MAP3K13(NP_001229246.1).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	108kDa
Uniprot ID	O43283
Dilution Range	WB: 1:1000-1:3000
Expiration Date	12 months from date of receipt.



Western blot
analysis of
extracts of
var...