



## Product Datasheet

### MAPK6 antibody (orb540332)


|                            |                                                         |
|----------------------------|---------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to MAPK6                     |
| <b>Reactivity</b>          | Human, Mouse, Rat                                       |
| <b>Conjugation</b>         | Unconjugated                                            |
| <b>Tested Applications</b> | ICC, IF, IHC, WB                                        |
| <b>Immunogen</b>           | Recombinant fusion protein of human MAPK6(NP_002739.1). |
| <b>Preservatives</b>       | PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.        |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.             |
| <b>Note</b>                | For research use only                                   |
| <b>Isotype</b>             | IgG                                                     |
| <b>Clonality</b>           | Polyclonal                                              |
| <b>MW</b>                  | 120kDa                                                  |
| <b>Uniprot ID</b>          | <a href="#">Q16659</a>                                  |
| <b>Dilution Range</b>      | WB: 1:200-1:2000                                        |
| <b>Expiration Date</b>     | 12 months from date of receipt.                         |



Immunofluorescence analysis of L929 cell...



Western blot analysis of extracts from n...



Western blot analysis of extracts of var...