


Product Datasheet

MAPK8IP1 antibody (orb540333)

Description	Rabbit polyclonal antibody to MAPK8IP1
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human MAPK8IP1(NP_005447.1).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	78kDa
Uniprot ID	Q9UQF2
Dilution Range	WB: 1:200-1:2000
Expiration Date	12 months from date of receipt.



Western blot
analysis of
extracts of
var...