

Product Datasheet

MS4A2 antibody (orb540477)

Description	Rabbit polyclonal antibody to MS4A2
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of human MS4A2(NP_000130.1).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	27kDa
Uniprot ID	Q01362
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.