



Product Datasheet

MTCO1 antibody (orb540488)

Description	Rabbit polyclonal antibody to MTCO1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	A synthetic peptide of human MTCO1(YP_002124304.1).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	40KDa/37KDa
Uniprot ID	P00395
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of NIH-3T3 C...



Western blot analysis of extracts of var...