



Product Datasheet

MTMR4 antibody (orb540500)


Description	Rabbit polyclonal antibody to MTMR4
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human MTMR4(NP_004678.3).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	133kDa
Uniprot ID	Q9NYA4
Dilution Range	WB: 1:200-1:2000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of U-2 OS cells...



Immunofluorescence analysis of L929 cells...



Western blot analysis of extracts of various tissues...