



Product Datasheet

Melatonin Receptor 1A antibody (orb540505)


Description	Rabbit Polyclonal antibody to Melatonin Receptor 1A
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human MTNR1A(NP_005949.1).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	39KDa
Uniprot ID	P48039
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:100-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Mouse br...



Immunohistochemical staining of Rat brai...



Western blot analysis of extracts of var...