



Product Datasheet

NDUFB7 antibody (orb540603)

Description	Rabbit polyclonal antibody to NDUFB7
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human NDUFB7(NP_004137.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	16kDa
Uniprot ID	P17568
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of L929 cell...



Western blot analysis of extracts of var...