



Product Datasheet

NMT1 antibody (orb540678)


Description	Rabbit polyclonal antibody to NMT1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human NMT1(NP_066565.1).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	57kDa
Uniprot ID	P30419
Dilution Range	WB: 1:200-1:2000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human ma...



Immunohistochemical staining of human co...



Immunohistochemical staining of rat ovar...