



Product Datasheet

NQO2 antibody (orb540715)


Description	Rabbit polyclonal antibody to NQO2
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human NQO2(NP_000895.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	26kDa
Uniprot ID	P16083
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Mouse kidney tissue.



Immunohistochemical staining of Human liver tissue.



Western blot analysis of extracts of various tissues.