



Product Datasheet

NUDT1 antibody (orb540756)


Description	Rabbit polyclonal antibody to NUDT1
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Recombinant fusion protein of human MTH1(NP_945192.1).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	18kDa
Uniprot ID	P36639
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200, IP: 1:20-1:50
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human gastric tissue.



Western blot analysis of extracts from 293T cells.



Western blot analysis of extracts from various cell lines.