


Product Datasheet

OTUD6B antibody (orb540813)

Description	Rabbit polyclonal antibody to OTUD6B
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human OTUD6B (NP_057107.3).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	38kDa
Uniprot ID	Q8N6M0
Entrez	51633
Dilution Range	WB: 1:500-1:2000, IP: 1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...