


Product Datasheet

PAFAH1B1 antibody (orb540836)

Description	Rabbit polyclonal antibody to PAFAH1B1
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human PAFAH1B1(NP_000421.1).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	47kDa
Uniprot ID	P43034
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot
analysis of
extracts of
var...