



Product Datasheet

PARP1 antibody (orb540853)


Description	Rabbit polyclonal antibody to PARP1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human PARP1(NP_001609.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	116kDa
Uniprot ID	P09874
Dilution Range	WB: 1:1000 - 1:2000, IHC-P: 1:50 - 1:200, IF/ICC: 1:50 - 1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...



Western blot analysis of extracts of var...



Immunofluorescence analysis of NIH/3T3 C...