



Product Datasheet

PCGF1 antibody (orb540873)


Description	Rabbit polyclonal antibody to PCGF1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human PCGF1(NP_116062.2).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	35KDa
Uniprot ID	Q9BSM1
Dilution Range	IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of U-2OS cel...



Immunofluorescence analysis of L-929 cel...



Western blot analysis of extracts of var...