


Product Datasheet

PGM2 antibody (orb540932)

Description	Rabbit polyclonal antibody to PGM2
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human PGM2(NP_060760.2).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	70kDa
Uniprot ID	Q96G03
Dilution Range	WB: 1:200-1:2000
Expiration Date	12 months from date of receipt.



Western blot
analysis of
extracts of
var...