


Product Datasheet

DPYSL2 (Phospho-T514) antibody (orb540951)

Description	Rabbit polyclonal antibody to DPYSL2 (Phospho-T514)
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around T514 of human CRMP2/DPYSL2(NP_001377.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	69kDa
Uniprot ID	Q16555
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...