



Product Datasheet

Histone H3 (Phospho-S28) antibody (orb540953)


Description	Rabbit polyclonal antibody to Histone H3 (Phospho-S28)
Reactivity	All, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	A synthetic phosphorylated peptide around S28 of human Histone H3(NP_003520.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	17kDa
Uniprot ID	Q16695
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of U2OS cell...



Immunofluorescence analysis of NIH/3T3 C...



Immunofluorescence analysis of H9C2 cell...