



Product Datasheet

IRE1alpha (Phospho-S724) antibody (orb540955)


Description	Rabbit polyclonal antibody to IRE1alpha (Phospho-S724)
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around S724 of human IRE1alpha(NP_001424.3?).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	126KDa
Uniprot ID	O75460
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of rat spleen...



Immunohistochemical staining of mouse spleen...



Western blot analysis of extracts of NIH3T3 cells.