

Product Datasheet

NPM1 (Phospho-T199) antibody (orb540960)

Description Rabbit polyclonal antibody to NPM1 (Phospho-T199)

Reactivity Human, Rat

Conjugation Unconjugated

Tested Applications WB

Immunogen A synthetic phosphorylated peptide around T199 of human NPM1(NP_001032827.1).

Preservatives PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.

Storage Store at -20°C. Avoid freeze / thaw cycles.

Note For research use only


Isotype IgG

MW 38kDa

Uniprot ID **P06748**

Dilution Range WB: 1:500-1:2000

Expiration Date 12 months from date of receipt.



Western blot analysis of extracts of var...