


Product Datasheet

P90RSK (Phospho-S363) antibody (orb540969)

Description	Rabbit polyclonal antibody to P90RSK (Phospho-S363)
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	A synthetic phosphorylated peptide around S363 of human p90Rsk/RSK1/RPS6KA1(NP_002944.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	90KDa
Uniprot ID	Q15418
Dilution Range	WB: 1:500-1:2000, IP: 1:50-1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...