


Product Datasheet

SRC (Phospho-Y529) antibody (orb540971)

Description	Rabbit polyclonal antibody to SRC (Phospho-Y529)
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	A synthetic phosphorylated peptide around Y529 of human Src(NP_005408.1).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	60KDa
Uniprot ID	P12931
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:100, IF: 1:100-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...