


## Product Datasheet

### PIK3C3 antibody (orb540982)

|                            |                                                                |
|----------------------------|----------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to PIK3C3                           |
| <b>Reactivity</b>          | Human, Mouse, Rat                                              |
| <b>Conjugation</b>         | Unconjugated                                                   |
| <b>Tested Applications</b> | ICC, IF, IHC, WB                                               |
| <b>Immunogen</b>           | Recombinant fusion protein of human PIK3C3/VPS34(NP_002638.2). |
| <b>Preservatives</b>       | PBS with 0.05% proclin300,50% glycerol,pH7.3.                  |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.                    |
| <b>Note</b>                | For research use only                                          |
| <b>Isotype</b>             | IgG                                                            |
| <b>Clonality</b>           | Polyclonal                                                     |
| <b>MW</b>                  | 100KDa                                                         |
| <b>Uniprot ID</b>          | <b>Q8NEB9</b>                                                  |
| <b>Dilution Range</b>      | WB: 1:500-1:2000                                               |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                |



Western blot analysis of extracts of var...