



Product Datasheet

PMS2 antibody (orb541035)

Description	Rabbit polyclonal antibody to PMS2
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IF, IP, WB
Immunogen	Recombinant fusion protein of human PMS2(NP_000526.2).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	110KDa
Uniprot ID	P54278
Dilution Range	WB: 1:500-1:2000, IF: 1:50-1:200, IP: 1:20-1:50
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of U2OS cell...



Western blot analysis of extracts of var...